

AFM Standards and References

Introduction

Standards and references are required for proper operation of an AFM and to verify optimal operating conditions and calibration of the instrument. Standards are used to assure that the absolute measurements are correct, while references assure that the instrument is giving consistent results.

For establishing calibrated topography measurements, the AFM scanner must be certified with calibration standards having pre-established dimensions in the X, Y and Z axis. References can be used for establishing that an AFM mode is operating correctly and for establishing the proper operation of the AFM.

This guide serves as a reference for finding standards and references that are helpful for making meaningful measurements with an AFM. When possible several vendors for a particular type of standard are listed.

Z Calibration

Calibration standards/references are needed to calibrate SPM in the vertical axis. For calibrations greater than 10 nm step height, standards or references are typically fabricated by etching patterns in a quartz substrate. Another source for references is etched silicon, or silicon dioxide coated with a uniform layer of metal. When calibrating the instrument for Z height measurements below 10 nm, nano-spheres, atomic terraces of silicon or HOPG may be used as a reference specimen.




FIGURE C-1 This AFM image is of a single etched pattern of an AFM standard. The histogram at the right show the height of the feature is 16.62 nm.

| Company | Standard / Reference |
|--|--|
| VLSI Standards www.vlsi.com | Standards Z=18nm, 44nm, 100nm, 180nm Pitch=1.8μm, 3μm, 5μm, 10μm, 20μm References: Z=18nm, 44nm, 100nm, 180nm Pitch=3μm, 10μm |
| MicroMasch www.spmtips.com | Z=20±25nm, 90±104nm, 496±503nm, 1000nm, 15000nm Pitch=3μm, 10μm HOPG all grades |
| SPI Supplies www.2spi.com | Pitch=300nm, 700nm Z~100nm |
| Ted Pella www.tedpella.com | Pitch=0.463μm Z~31nm |
| Pacific Nanotechnology www.pacificnanotech.com | Pitch=20μm, 10μm, 5μm, 3μm Z~75nm |
| Electron Microscopy Science www.emsdiasum.com | Pitch=150nm, 300nm, 700nm Z~100nm |
| NTT-AT www.ntt-at.com | AtomicStep=3.1nm |

XY Linearity / Calibration / Orthogonality

Calibration of the AFM in the XY axis so that meaningful dimension measurements can be made requires establishing the linearity and orthogonality of the AFM scanner. This is typically done with etched

Appendix C • Standards and References

patterns in quartz or silicon. Often the specimen that is used for calibrating the Z axis can be used for the XY axis as well.

| Company | Standard / Reference |
|--|---|
| VLSI Standards www.vlsi.com | Standards Z=18nm, 44nm, 100nm, 180nm Pitch=1.8µm, 3µm, 5µm ,10µm, 20µm References: Z=18nm, 44nm, 100nm, 180nm Pitch=3µm,10µm |
| MicroMasch www.spmtips.com | Pitch=3µm, 10µm Z~1µm, 2µm |
| Pacific Nanotechnology www.pacificnanotech.com | Pitch=20µm, 10µm, 5µm, 3µm Z~75nm |

XZ, YZ Orthogonality

An AFM may be calibrated in XY and calibrated in Z and may not be useful for making angular measurements. This is because the XZ and YZ axis may not be orthogonal. With the orthogonality references, this problem can be avoided. This reference is fabricated by making a 1-D array, or line, of triangles in a silicon wafer.

| Company | Standard / Reference |
|--|---|
| MicroMasch www.spmtips.com | Pitch=3.0µm ±5nm Triangle edge curvature<10nm Top angle=70deg |




FIGURE C-2 This AFM image of a triangle pattern that is etched in silicon is used for measuring the cross talk in the XZ and YZ axis of the scanner. Measurements from the line profile show that the angles are 55.16 and 55.59 degrees on the left and the right side of one of the triangular features.

Flatness / Roughness


An AFM is capable of measuring surface roughness of samples at the nanometer scale. However, it is important to establish that the surface roughness measured with the AFM is not the noise floor of the instrument. Also, when measuring flat samples, it is important to know the “bow” associated with the scanner. These references are useful for establishing the performance of the instrument with respect to noise floor and bow.

| Company | Standard / Reference |
|--|-----------------------------|
| NT-MDT www.ntmdt.com | Sa=0.06nm |
| SPI Supplies www.2spi.com | Atomic corrugations on HOPG |

Particles

SPM height calibration can be done on nanospheres of uniform size; it also verifies SPM performance on samples with weak adhesion to the surface. Often it is helpful to include nanoparticles with a sample to establish the sizes of features in the imager relative to the nanoparticles.


| Company | Standard / Reference |
|--|---|
| Duke Scientific Co. www.dukescientific.com | 20<D<900nm Traceable and certified polymer spheres come in dry powder and aqueous suspensions in a great variety of sizes. |
| TedPella Inc. www.tedpella.com | D=5nm D=15nm D=30nm AFM gold standard kit consists of three sizes of gold colloids in aqueous suspensions. |
| American Dye Source Inc. www.adsdyes.com EvidentTechnology Inc. www.evidenttech.com | All colors of q-dots. Quantum dots of any size come in great variety of kits |



Bio-Molecules

Often it is necessary to verify SPM image quality on soft samples or samples with weak adhesion to the surface. It is important that this type of reference have a long-life, be stable and be indestructible. Before measuring images of “unknown” biological molecules it is often helpful to practice on this “known” sample.

| Company | Standard / Reference |
|--|------------------------|
| NTI-Europe www.ntinstruments.com | Z=2nm Length~1009nm |



Friction Mode

LFM reference comes as 1-D array of triangular steps having precise linear and angular dimensions. Establishing quantitative LFM data requires standardizing the AFM scanner output in the vertical and horizontal axis.

| Company | Standard / Reference |
|--|---|
| MicroMasch www.spmtips.com | Pitch=3.0 μ m \pm 5nm Edge curvature<10nm Top angle=70deg |

Vibrating Phase Mode


In vibrating mode, the AFM can measure changes in chemical composition of polymer samples. A reference of polymer material that has hard and soft regions is used for this purpose.

| Company | Standard / Reference |
|---|----------------------|
| Electron Microscopy Science www.emsdiasum.com | As small as 10nm |

Tip Visualization

Tip Visualization helps visualize the geometry of the scanning probe without the costly need of SEM. Tip radius sample exhibits features sharper than an AFM tip. First-order approximation of the tip apex can be obtained.

| Company | Standard / Reference |
|--|----------------------|
| Aurora NanoDevice Inc. www.aurorand.com | Z~30nm X,Y~10nm |



Conclusion

Making consistent, or even quantitative, measurements with an AFM requires diligence and the appropriate standard or reference. This partial list of standards/references can be used for making consistent or quantitative measurements.